
	<b>GILES CHEMICAL ~ PREMIER MAGNESIA</b>		
	<b>Company Procedure</b>		
	Title: <b>Steps for Liquid Load Testing</b>	Number: <b>L12-PR-100-024</b>	
	Owner: <b>Ashley Williams</b>	Revision: <b>3</b>	
	Effective Date: <b>04/29/16</b>	Page: <b>1 of 3</b>	

## 1.0 Purpose

The purpose of this procedure is to define the steps for testing liquid loads.

## 2.0 Scope

This procedure applies to all finished liquid loads.

## 3.0 Responsibility

QA Personnel are responsible for liquid load testing.

## 4.0 Safety Considerations

Proper PPE is required while in the lab, including but not limited to steel-toed shoes and safety glasses. Safety is a condition of employment. Employees are not authorized to work in an unsafe manner and are prohibited from harming the environment of the facility or community.

## 5.0 Materials/Equipment

- 500mL Beaker
- Magnetic Stirring Bar – 1” length
- Microwave
- De-ionized Water
- pH Meter System
- Magnetic Stirring Plate
- Weighing Balance (accurate to 0.01 or higher)
- 25mL Flask
- Aqua Tester Comparator
- 2 Aqua Tester Comparator Tubes with glass caps
- Calculator
- Thermometer



## 6.0 Procedure

All liquid loads are to be tested for specific gravity, pH, and color.

Liquid Loads must be heated, using the microwave, until all crystallized salt is fully dissolved in the solution.

### Controlled Document

Only those quality documents viewed through the Giles Chemical electronic Documentation System are officially controlled. All other copies, whether viewed through another computer program or a printed version, are not controlled and, therefore, the Quality Unit at Giles assumes no responsibility for accuracy of the document.

	GILES CHEMICAL ~ PREMIER MAGNESIA		
	Company Procedure		
	Title: Steps for Liquid Load Testing	Number: L12-PR-100-024	
	Owner: Ashley Williams	Revision: 3	
	Effective Date: 04/29/16	Page: 2 of 3	

## A. pH - Temperature

### Procedure:

1. Place liquid sample in 500ml beaker with stir bar.
  2. Heat sample to until all crystallized salt is in solution.
  3. Remove the protective cap from the probe tip
  4. Rinse the tip with de-ionized water and wipe with Chem-wipe or paper towel.
  5. Place the probe in the sample and turn on the pH meter.
  6. Turn on the magnetic stirring plate and stir the sample to ensure an accurate reading.
  7. Record the pH value and temperature when pH meter reads ready.
  8. Record pH and temperature on *Final Product Liquid – Daily Quality Control Report (L12-FM-100-006)*.
  9. Remove probe, rinse with de-ionized water, dry the probe and replace the protective cap.
- **The instrument can be turned on and off with out losing its calibration. The meter should be calibrated daily and anytime the meter shows an obviously erroneous reading.**

## B. Specific Gravity-

### Procedure:



1. A dry 25 ml volumetric flask is placed on the weighing balance and tarred to zero.
2. 25 ml of subject sample is added to the volumetric flask and the weight recorded.
3. Specific gravity is determined using the following formula

$$\frac{\text{Weight of sample (g)}}{\text{Volume of sample (mL)}} = \text{Specific Gravity (g / mL)}$$

4. % MgSO<sub>4</sub> is determined by using the Percent MgSO<sub>4</sub> chart.(Chart is located in the cabinet drawer below the balance).
5. Record result on *Final Product Liquid – Daily Quality Control Report (L12-FM-100-006)*.

### Controlled Document

Only those quality documents viewed through the Giles Chemical electronic Documentation System are officially controlled. All other copies, whether viewed through another computer program or a printed version, are not controlled and, therefore, the Quality Unit at Giles assumes no responsibility for accuracy of the document.

	GILES CHEMICAL ~ PREMIER MAGNESIA		
	Company Procedure		
	Title: Steps for Liquid Load Testing	Number: L12-PR-100-024	
	Owner: Ashley Williams	Revision: 3	
	Effective Date: 04/29/16	Page: 3 of 3	

## C. Color Test

### Procedure:

1. Fill one of the tubes in the aqua tester comparator to the mark with sample.
2. Place the glass cap piece on the top of the tube.
3. Place the tube in the right hand side of the comparator tube holder.
4. Fill the other tube to the mark with de-ionized water.
5. Place the glass cap piece on the top of the tube.
6. Place the tube in the left hand side of the comparator tube holder.
7. Turn color dial until both tubes look the same color in the comparator.
8. Record result on liquid *Final Product Liquid – Daily Quality Control Report (L12-FM-100-006)*.

## 7.0 Reference Documents

*Final Product Liquid – Daily Quality Control Report (L12-FM-100-006)*

## 8.0 Change Information

Updated to new SOP template. Added use of thermometer.

### Controlled Document

Only those quality documents viewed through the Giles Chemical electronic Documentation System are officially controlled. All other copies, whether viewed through another computer program or a printed version, are not controlled and, therefore, the Quality Unit at Giles assumes no responsibility for accuracy of the document.