

Safety: Wear the appropriate PPE when working in the Lab.

Purpose: Determination of Bulk Density for Crystalline Salt Samples

Procedure:

Background Information:

It is sometimes important to know how much space a certain weight of a crystalline product will take up in transport and/or storage. This is where the term bulk density is used. Bulk density for crystalline samples is analogous to the density for liquid samples. It denotes the weight of a certain volume of crystalline sample or the volume of a certain weight of a crystalline sample. Bulk density is an intrinsic property of the material being tested and does not depend on the amount of the sample present. The following procedure is used to determine bulk density.

Scope:

A crystalline sample is thoroughly compacted to a known volume in a tarred graduated cylinder. The mass of the sample in the cylinder is measured. The bulk density is calculated and converted to lb/ft³.

Equipment:

100-mL graduated cylinder Funnel Small spatula Weighing Balance

- 1. A 100-mL graduated cylinder is placed on the balance and tared to zero.
- 2. The Cylinder is filled with the salt sample until the level reaches 100ml.
- 3. The cylinder is tapped lightly with spatula until sample ceases to drop.
- 4. The volume of the sample in the cylinder is recorded.
- 5. The cylinder is placed back on the balance and weight of the sample is recorded.
- 6. The bulk density can be determined by using the following calculation

weight of salt sample (g) / Volume of salt sample (mL) = bulk density (g/mL) g / mL x 62.4 = lb. / ft³



GILES CHEMICAL

Determination of Bulk Density for Crystalline Salt Page . 2

Page : 2 of 3

Revision Date

01 06/29/2009

Author: Carl Mooney Job Specific

TRAINING DOCUMENTATION

| | EMPLOYEE | TITLE | SIGNATURE | DATE |
|----|----------|-------|-----------|------|
| 1 | | | | |
| 2 | | | | |
| 3 | | | | |
| 4 | | | | |
| 5 | | | | |
| 6 | | | | |
| 7 | | | | |
| 8 | | | | |
| 9 | | | | |
| 10 | | | | |
| 11 | | | | |
| 12 | | | | |
| 13 | | | | |
| 14 | | | | |
| 15 | | | | |
| 16 | | | | |
| 17 | | | | |
| 18 | | | | |
| 19 | | | | |
| 20 | | | | |
| 21 | | | | |
| 22 | | | | |
| 23 | | | | |
| 24 | | | | |
| 25 | | | | |
| 26 | | | | |
| 27 | | | | |
| 28 | | | | |
| 29 | | | | |
| 30 | | | | |



GILES CHEMICAL

COMPANY PROCEDURE

Determination of Bulk Density for Crystalline Salt | Page | 2

Page : 3 of 3

Revision Date 01 06/29/2009

Author: Carl Mooney Job Specific

| Revision Number | Revision Date | Revision Author | Revision Description |
|--------------------|------------------|--------------------|---|
| 00 | 03/28/06 | СМ | New Document |
| 01 | 06/29/2009 | SL | -Remove compacting the sample initially -3 year revision -Placed on new form. |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |