
	GILES CHEMICAL ~ PREMIER MAGNESIA		
	Company Procedure		
	Title: Analysis of Filter Press Cake	Number: L12-PR-100-032	
	Owner: Hunter Douglas	Revision: 2	
	Effective Date: 05/15/2015	Page: 1 of 5	

1.0 Purpose:

The purpose of this procedure is to provide a method for testing the filter press cake in order to determine efficiency and hint at any deviations to the normal process.

2.0 Scope:

This procedure is to be performed at intervals specified by management and when a deviation to the normal process occurs.

3.0 Responsibility:

QA personnel or representative will request a sample from the material handler and will provide the material handler with a sample bag that has a blank label with spots for, at minimum, the date and the MgO railcar #.

This procedure is to be performed by any qualified laboratory personnel; a second analyst will review data for accuracy and completeness.

4.0 Safety Considerations:

Proper PPE should be worn at all times during this procedure. Including but not limited to gloves, safety goggles, and lab coat.

Safety is a condition of employment. Employees are not authorized to work in an unsafe manner and are prohibited from harming the environment of the facility or community.



5.0 Materials/Equipment:

Equipment:

- Large crucibles
- Drying oven (100 ° C)
- Forceps
- Desiccator
- Vacuum filter assembly including Buchner funnel, stopper, tubing and filter pump
- Whatman grade 40 (8 µm) filter paper – OR – other low retention filter paper
- Balance (0.0000g accuracy min.)
- Weigh Paper
- Spatula
- pH Meter
- 1000-ml Volumetric Flask

Controlled Document

Only those quality documents viewed through the Giles Chemical electronic Documentation System are officially controlled. All other copies, whether viewed through another computer program or a printed version, are not controlled and, therefore, the Quality Unit at Giles assumes no responsibility for the accuracy of the document.

	GILES CHEMICAL ~ PREMIER MAGNESIA		
	Company Procedure		
	Title: Analysis of Filter Press Cake	Number: L12-PR-100-032	
	Owner: Hunter Douglas	Revision: 2	
	Effective Date: 05/15/2015	Page: 2 of 5	

- 100-mL Volumetric Flask
- 100-mL Graduated Cylinder
- 250-mL Filter Flasks
- Watch Glass
- Stir bar
- Stir plate
- Class A – 100 ml burette
- Burette Stand
- 1000-μL Eppendorf Pipette and Tips
- 5-mL Eppendorf Pipette and Tips

Reagents:

- Ammonium Chloride
- DI H₂O
- Eriochrome Black TS
- 0.05 M Edetate Disodium (Disodium EDTA) Volumetric Solution
- 3 N Hydrochloric Acid Solution (If Needed)
- 1 N Sodium Hydroxide Solution (If Needed)
- Lithium Metaborate
- Lithium Bromine
- Nitric Acid (conc.)

6.0 Procedure:

Solutions Preparation: (If Needed)

Ammonium – Ammonium Chloride Buffer Test Solution– Dissolve 67.5 g of ammonium chloride in water, add 570 ml of ammonium hydroxide, dilute with DI H₂O to 1000 ml.

2% Nitric Acid (HNO₃) – Add 57.1mL of conc HNO₃ to a 2000mL volumetric flask that is half filled with DiH₂O. Dilute to volume with DiH₂O and mix well by inversion.

10% Nitric Acid (HNO₃) – Add 143mL of conc HNO₃ to a 1000mL volumetric flask that is half filled with DiH₂O. Dilute to volume with DiH₂O and mix well by inversion.



Sample Prep/Analysis:

LOD:

1. From the press cake sample received, cut a lab sample perpendicular to the cake surface weighing about 30 grams and proceed to analyze each individually as directed below.

Controlled Document

Only those quality documents viewed through the Giles Chemical electronic Documentation System are officially controlled. All other copies, whether viewed through another computer program or a printed version, are not controlled and, therefore, the Quality Unit at Giles assumes no responsibility for the accuracy of the document.

	GILES CHEMICAL ~ PREMIER MAGNESIA		
	Company Procedure		
	Title: Analysis of Filter Press Cake	Number: L12-PR-100-032	
	Owner: Hunter Douglas	Revision: 2	
	Effective Date: 05/15/2015	Page: 3 of 5	

2. Weigh a 100 mL crucible. **Record the weight (A).**
3. Place the lab sample in the crucible, weigh, and **record the weight (B)**

B – A = the weight of the sample

4. Place the crucible and sample in the 100°C oven overnight.
5. Next day, using the forceps, remove the crucible from the oven and place in the desiccator.
6. When cool, weigh the crucible containing the sample, and **record the weight (D).**

D – A ÷ (B-A) · 100% = % Total Solids in the Filter Press Cake

B – D ÷ (B-A) · 100% = % LOD of the Filter Press Cake

MgSO₄ Assay:

1. Grind the entire dry portion and weigh out a 0.50 g portion and **record weight (C)**. Then place the portion into a suitable beaker.
2. Add exactly 50mL of DI water to the beaker, add a stir bar, and stir until dissolved.
3. Using a Whatman grade 40 (8 µm) filter paper, filter the sample through a vacuum and rinse the beaker, stir bar, and filtrate with another 50mL portion of DI water.
5. Titrate the filtrate using the following procedure *USP Assay: Magnesium Sulfate (L12-PR-100-008) (steps 7 to 14)*.
6. Use the formula below to calculate % MgSO₄ in the Filter Press Cake:

$$\frac{\text{mL of EDTA solution} \cdot 0.006018}{C} \cdot 100\% = \% \text{ MgSO}_4 \text{ in the Filter Press Cake (D)}$$



Quantification of Unreacted Matter:

1. Accurately weigh ~50mg of the dry grinds out onto a sheet of weigh paper.
2. Using the KATANA fluxer, flux the 50mg with 1.0g of **Lithium Metaborate**, and 250mg of **Lithium Bromine** into 100mL of 10% HNO₃.
(Note: Ensure the sample is entirely in solution before continuing.)
3. Prepare a series of ICP standards in the following concentrations:

Standard ID	Ca Conc.	Fe Conc.	Mg Conc.
STD 1	0	0 ppm	0 ppm
STD 2	2 ppm	1 ppm	5 ppm
STD 3	4 ppm	2 ppm	12.5 ppm
STD 4	6 ppm	3 ppm	20 ppm

Controlled Document

Only those quality documents viewed through the Giles Chemical electronic Documentation System are officially controlled. All other copies, whether viewed through another computer program or a printed version, are not controlled and, therefore, the Quality Unit at Giles assumes no responsibility for the accuracy of the document.



	GILES CHEMICAL ~ PREMIER MAGNESIA		
	Company Procedure		
	Title: Analysis of Filter Press Cake	Number: L12-PR-100-032	
	Owner: Hunter Douglas	Revision: 2	
	Effective Date: 05/15/2015	Page: 4 of 5	

4. Using the method labeled Mud Analysis in the SALSA ICP software, initiate the ICP as directed in *USP ICP-OES Analysis (L13-PR-100-057)*.
5. Place standards in positions 8-11. Align the torch using the *Auto Align* function and STD 4. Repeat the *Auto Align* function until the window shows a good saturation and a dx/dy of 0/0.
6. Prepare the sample by diluting the fluxed solution by 1:10. (*Example: Pipette 10mL into a 100mL volumetric flask and q.s. with 2% HNO₃*)
7. Under the **Method**→**QC Automation** tab ensure that standards 1-4 and QC1 are selected. Continue to the **Analysis** tab and create a new analysis.
8. Continue to the **Sequence** pane and start the sequence.
9. After the standards have run check the calibration curve under **Methods**→**Elemental Selection**→**Calibration**. (*Note: The line for each element should be relatively linear and should have a Rho value near 1.0*)
10. If the calibration curve and the QC pass then proceed otherwise use the **Auto Align** function again and re-establish a new calibration curve. (*Note: To determine if the QC passes go to the Analysis*→*Results* *tab. If beside the concentrations of the elements there is an (L) or (H) then the QC failed and must be re-run, if there is no such character then continue.*)
11. Go back to the **Method**→**QC Automation** tab and uncheck the standards and the first QC. Continue to the **Sequence** pane and check the number of boxes = to the number of samples to be analyzed. Give each sample a name that includes the date the sample was taken. In the weight cell enter 0.050 and in the volume cell enter 1000. (*Note: The values entered are the values from the sample prep i.e.- {[50mg/100mL] x [10mL/100mL]} if the sample was prepped differently then the numbers will change.*)
12. Click the **Update Sequence** button and then **Run Sequence**. There should be a QC standard run every 5 samples and at the end of the run.
13. Once the run is complete go to the **Analysis** pane again and go to the **Repots** tab. Select all of the injections from the first QC down. Then click the **Load Report** button. Select **Simple SBD** report and the **Format 2** button. Click **Print** and the **Printer** button. Ensure the printer has paper loaded.
14. The report will show values in the concentration cell. These values are reported in ppm so to convert to % the value will need to be divided by 10,000.
15. To get % Mg⁺² the Mg associated with the MgSO₄ must be subtracted using the calculation:

$$\frac{Mg \text{ (ppm)}}{10000} - (D * 0.2020) = \% Mg^{+2}$$

Controlled Document

Only those quality documents viewed through the Giles Chemical electronic Documentation System are officially controlled. All other copies, whether viewed through another computer program or a printed version, are not controlled and, therefore, the Quality Unit at Giles assumes no responsibility for the accuracy of the document.

	GILES CHEMICAL ~ PREMIER MAGNESIA		
	Company Procedure		
	Title: Analysis of Filter Press Cake	Number: L12-PR-100-032	
	Owner: Hunter Douglas	Revision: 2	
	Effective Date: 05/15/2015	Page: 5 of 5	

7.0 Reference Documents:

USP Assay: Magnesium Sulfate (L12-PR-100-008)

USP ICP-OES Analysis (L13-PR-100-057)

Analysis of Filter Cake (L15-PR-100-F064)

8.0 Change Information:

Update of entire document

Controlled Document

Only those quality documents viewed through the Giles Chemical electronic Documentation System are officially controlled. All other copies, whether viewed through another computer program or a printed version, are not controlled and, therefore, the Quality Unit at Giles assumes no responsibility for the accuracy of the document.