



|  |                                   |                               |   |
|--|-----------------------------------|-------------------------------|---|
|  | GILES CHEMICAL ~ PREMIER MAGNESIA |                               |  |
|  | Company Procedure                 |                               |   |
|  | Title: <b>USP pH</b>              | Number: <b>L12-PR-100-002</b> |   |
|  | Owner: <b>Stephen Ballew</b>      | Revision: <b>1</b>            |   |
|  | Effective Date: <b>03/31/13</b>   | Page: <b>1 of 3</b>           |   |

## 1.0 Purpose

To describe how to determine the pH of magnesium sulfate heptahydrate in solution following USP Monograph: Magnesium Sulfate, and General Chapter <791>.

## 2.0 Scope

This procedure applies to USP lot change, stability testing, and any time USP quality needs to be verified. All USP testing is performed in the Quality Assurance laboratory.

## 3.0 Responsibility

QA Lab personnel are responsible for USP testing.

## 4.0 Safety Considerations

Wear safety glasses and lab coat when working in the lab.

Safety is a condition of employment. Employees are not authorized to work in an unsafe manner and are prohibited from harming the environment of the facility or community.

## 5.0 Materials/Equipment

- pH Meter System – VWR SB-20
- Balance – Mettler Toledo X5105Du, B13929Z316
- 150 ml Beaker
- Supply of De-ionized Water
- Supply of Buffer Solutions – pH 4.00, 7.00 and 10.00
- Magnetic Stirring Plate
- Magnetic Stirring Bars – 1” length
- Sheet of 8½” x 11” office letter paper



## 6.0 Procedure

**The pH system should be calibrated once per working day.**

1. (NOTE: If using *Test Solution* from USP Identification Test, skip to step 3). Weigh approximately 5.00 g of the sample on a piece of paper using the analytical balance.
2. To a 150 ml beaker add 100 ml of H<sub>2</sub>O.
3. Place the beaker on the magnetic stirring plate and add the sample.
4. Add a stir bar to the beaker and turn stirring plate on (½ to ¾ max. setting).

### Controlled Document

Only those quality documents viewed through the Giles Chemical electronic Documentation System are officially controlled. All other copies, whether viewed through another computer program or a printed version, are not controlled and, therefore, the Quality Unit at Giles assumes no responsibility for accuracy of the document

|  |                                   |                               |   |
|--|-----------------------------------|-------------------------------|---|
|  | GILES CHEMICAL ~ PREMIER MAGNESIA |                               |  |
|  | Company Procedure                 |                               |   |
|  | Title: <b>USP pH</b>              | Number: <b>L12-PR-100-002</b> |   |
|  | Owner: <b>Stephen Ballew</b>      | Revision: <b>1</b>            |   |
|  | Effective Date: <b>03/31/13</b>   | Page: <b>2 of 3</b>           |   |

- Remove protective laboratory film from the probe tip, remove probe tip from storage solution, rinse the tip with de-ionized water, and wipe with paper towel.
- When the sample has completely dissolved, place the probe in the solution and turn on the pH meter.
- Record the pH value one minute later. The level will have stabilized sufficiently in that interval.
- Remove probe, rinse with de-ionized water, place back into storage solution and replace the protective laboratory film.

**The pH of magnesium sulfate heptahydrate should be between 5.0 and 9.2 in a 1 to 20 solution.**

### **Calibration and Maintenance of the VWR Model SB20 pH meter (daily and weekly):**

#### **Daily Calibration:**

- Place the pH electrode into the pH 4.0 calibration buffer.
- Press lower left key on the pH meter to turn the unit on.
- Press the cal key and "calibrate" will be displayed in the lower field. "P1" will be displayed indicating it is ready for the first buffer point.
- When "ready" flashes, record that value on pH Meter Calibration Log. And press "ok" to accept the pH value. "P2" will then be displayed.
- Rinse the electrode with distilled water and place the electrode into the pH 7.0 buffer.
- When "ready" flashes again, record that value on pH Meter Calibration Log and press the "ok" key. "P3" will then be displayed.
- Rinse electrode with distilled water and place the electrode into the pH 10.0 buffer.
- When "ready" flashes again, record that value on pH Meter Calibration Log and press the "ok" key.
- The display will freeze for 2 seconds, and then the slope will momentarily be displayed. Record the slope on pH Meter Calibration Log. The slope should be between 92-102%.
- Meter will automatically advance to the measure mode. Rinse electrode and place into sample or storage solution.



#### **Weekly Maintenance:**

##### **To be done before daily calibration**

- Inspect the electrode for scratches, cracks, salt/crystal build up, or membrane/junction deposits.
- Rinse off salt build up with distilled water and remove any other deposits as indicated in manual.
- Replace storage solution.
- Record maintenance completion on pH Meter Calibration Log along with calibration data.

#### **Controlled Document**

Only those quality documents viewed through the Giles Chemical electronic Documentation System are officially controlled. All other copies, whether viewed through another computer program or a printed version, are not controlled and, therefore, the Quality Unit at Giles assumes no responsibility for accuracy of the document

|   |                                   |                        |   |
|---|-----------------------------------|------------------------|---|
|  | GILES CHEMICAL ~ PREMIER MAGNESIA |                        |  |
|   | Company Procedure                 |                        |   |
|   | Title: USP pH                     | Number: L12-PR-100-002 |   |
|   | Owner: Stephen Ballew             | Revision: 1            |   |
|   | Effective Date: 03/31/13          | Page: 3 of 3           |   |

## 7.0 Reference Documents

N/A

## 8.0 Change Information

Updated using *SOP Template Instructions (Q12-PR-100-004)* and *Document Numbering (Q12-PR-100-003)*

### Controlled Document

Only those quality documents viewed through the Giles Chemical electronic Documentation System are officially controlled. All other copies, whether viewed through another computer program or a printed version, are not controlled and, therefore, the Quality Unit at Giles assumes no responsibility for accuracy of the document